

Urinary Catheter Securement Device

PRODUCT DATA SHEET



DESCRIPTION

The urinary catheter securement device provides securement of the catheter in two simple steps; adhesive strip ensures the catheter does not slide once in position and the soft material tabs cross over each other to secure the catheter in place.

Available in two sizes; standard which will suit all 2-Way Urinary Catheters and a larger size to suit 3-Way Urinary Catheters.

FEATURES AND BENEFITS

- Two stabilisation steps ensure maximum securement to prevent sliding and dislodgement
- Low irritant, soft adhesive dressing component
- Soft material tabs can be re-applied numerous times to reposition
- Easy to apply and remove
- Standard size is suitable for 2-Way Urinary Catheters sized 8-26Fr
- Large size is suitable for 3-Way Urinary Catheters
- Upon removal, skin remains free from residue and does not require any additional products to assist with removal

PRODUCT DATA

- Standard size: 6.5cm x 10.5 cm
- Large size: 8.5cm x 12.5cm
- MRI compatible: Yes
- Latex free
- Materials: Velcro, non-woven fabric, acrylic adhesive, glassine paper
- Shelf Life: 5 Years
- ARTG number: 299241

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
UR060000	Standard	50 pcs
UR060001	Large	50 pcs